



Transperineal Needle Guide (GTK155)

Contact:

sales@usmedicalinn.com

www.usmedicalinn.com

+1 (305) 772-3806





GTK155

Transperineal biopsy needle guide Experience precision and safety with the GTK155 Transperineal Biopsy Needle Guide. Designed specifically for use with BK 8848 and 9048

Features:

- **Minimally Invasive Freehand Approach:** Enables a freehand, ultrasound-guided technique for transperineal prostate biopsies, offering flexibility and precision in needle placement.
- **Enhanced Patient Safety:** The disposable design minimizes the risk of infection and cross-contamination, ensuring a sterile procedure every time. This approach can reduce the need for antibiotics, lowering the risk of antibiotic resistance.
- **Comprehensive Prostate Access:** Compatible with 15G-18G instruments, the guide allows improved access to all prostate zones, including difficult-to-reach anterior areas, increasing the detection rate of clinically significant diseases.

GTK155

Compatible probes



BK:

BK8848

BK9048





GTK155

Ordering information

US Medical Catalog Number	Description
GTK155	Needle Guide Only



sales@usmedicalinn.com

www.usmedicalinn.com

+1 (305) 772-3806